



## DANAGENE SPIN SALIVA DNA KIT

Ref.0603.SPIN50 50 PREPS  
 Ref.0603.SPIN250 250 PREPS

### 1. INTRODUCTION

**DANAGENE SPIN SALIVA DNA kit** is designed for the rapid purification of **highly pure genomic DNA from saliva samples** using silica-based system with a MicroSpin format:

a) Preserved saliva samples in the **DANASALIVA Sample Collection Kit**.



b) Fresh saliva samples.

### 2. KIT COMPONENTS

|                        | 50 PREPS  | 250 PREPS | T <sup>a</sup> Stock |
|------------------------|-----------|-----------|----------------------|
| Lysis Buffer PS        | 35 ml     | 160 ml    | Room temperature     |
| Desinhibition Buffer * | 18 ml     | 82.50 ml  | Room temperature     |
| Wash Buffer *          | 10 ml     | 50 ml     | Room temperature     |
| Elution Bufer          | 6 ml      | 30 ml     | Room temperature     |
| Proteinase K*          | 30 mg     | 5 x 30 mg | -20°C                |
| MicroSpin columnas     | 50 units  | 250 units | Room temperature     |
| Collection Tubes       | 100 units | 500 units | Room temperature     |

(\* ) These solution must be prepared as indicated in the Preliminary Preparations section of the protocol.

*PRECAUTIONS: The Desinhibition Buffer contains guanidium salts, and should be handled with care. Guanidinium salts form highly reactive compounds when combined bleach.*

#### Intended Use

All DANAGENE products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

## 2.1 Equipment and additional reagents required

### Reagents

- 96 – 100 % ethanol

### Consumables

- 1.5 mL microcentrifuge tubes
- Disposable pipette tips

### Equipment

- Manual pipettors
- Heat block, dry bath, or water bath (70°C)
- Centrifuge for microcentrifuge tubes
- Personal protection equipment (e.g., lab coat, gloves, goggles)

## 2.1 Storage and stability

All components are stable for 12 months from the date of purchase being stored correctly and at room temperature (15-25°C).

## 3. PROTOCOL

### 3.1 Preliminary Preparations

- Dissolve the proteinase K in **1.3 ml of nuclease-free water** and store at –20°C. It is recommended to do several aliquots to avoid many thaw/freeze cycles. At this temperature it is stable for 1 year.
- **Add 10 ml (50 preps) or 50 ml (250 preps) of Ethanol 100 %** to the Desinhibition Buffer. Keep the container closed to avoid the ethanol evaporation.
- **Add 40 ml (50 preps) or 200 ml (250 preps) of Ethanol 100 %** to the Wash Buffer. Keep the container closed to avoid the ethanol evaporation.
- Pre-heat the Elution Buffer at 70°C.

### 3.2 Isolation from preserved saliva samples with the DANASALIVA Sample Collection Kit ( Ref. 0603.43)

1. Vortex DANASALIVA Sample Collection Kit tube containing preserved saliva sample in order to homogenate the sample correctly. Transfer **400 µl in a 1.5 ml microtube**.
2. **Add 600 µl Lysis Buffer PS + 25µl Proteinase K**. Vortex.
3. Incubate **at 55°C for 20 minutes** with vortex periodically.
4. Centrifuge **at 14.000 rpm for 2 minutes**.
5. Transfer **supernatant** in a new 1.5 ml microtube avoiding touching the pellet.
6. **Add 300 µl of Ethanol 100%**. Mix well.
7. Load **650 µl** mixture sample into reservoir of a combined MicroSpin column–collection tube assembly. **Centrifuge at 8.000 rpm for 30 seconds**.

8. Discard the flow-through and place the Spin Column back into the same 2 ml Collection Tube, repeat step 6 with the remaining sample mixture.
  9. Carefully place the Spin Column into a clean 2 ml Collection Tube (provided). Avoid splashing any flow-through onto the Spin Column.
  10. Add **500 µl of Desinhibition Buffer. Centrifuge at 12.000 rpm for 1 minute.** Discard the flow-through.
  11. Add **700 µl of Wash Buffer. Centrifuge at 14.000 rpm for 1 minute.** Discard the flow-through.
  12. **Dry silica membrane.** Centrifuge at 14.000 rpm for 3 minutes.
12. Place the MicroSpin Column into a 1.5 mL nuclease-free tube (not provided) and add **50µL Pre-heat the Elution Buffer** at 70°C on the centre of the white filter membrane. Incubate **at room temperature for 2 minutes.**
13. **Centrifuge** the spin column-tube assembly **at 14.000 rpm for 1 minute**, then discard the column. The DNA is now ready for downstream applications.

### **3.3 Protocol for genomic DNA isolation from 600-800 µl of fresh saliva**

**NOTE:** Do **NOT** eat, drink, smoke or chew gum for 30 minutes before giving your saliva sample.

Fresh saliva samples must be processed immediately or keep at 4°C if they will be processed in less than 2 hours.

1. Centrifuge **600-800 µl of saliva at 14.000 rpm** for 90 seconds. Remove the supernatant with micropipette without damaging the visible white pellet of cells.

If the cell pellet is very small, you can add another **600 µl of saliva** and repeat point 1.

Genomic DNA yield is proportional to the size of the cell pellet, the larger the pellets, the greater the amount of DNA obtained.

2. **Add 600 µl Lysis Buffer PS + 25µl Proteinase K.** Resuspend the cell pellet well using a micropipette. **Incubate at 55°C for 20 minutes** with vortex periodically.

3. Centrifuge **at 14.000 rpm for 2 minutes.**

4. Transfer **supernatant** in a new 1.5 ml microtube avoiding touching the pellet.

**5. Add 100 µl of Ethanol 100%.** Mix well.

6. Load lysate mixture sample into reservoir of a combined MicroSpin column-collection tube assembly. **Centrifuge at 8.000 rpm for 30 seconds.**

7. Carefully place the Spin Column into a clean 2 ml Collection Tube (provided). Avoid splashing any flow-through onto the Spin Column.
8. Add **500 µl of Desinhibition Buffer. Centrifuge at 12.000 rpm for 1 minute.** Discard the flow-through.
9. Add **700 µl of Wash Buffer. Centrifuge at 14.000 rpm for 1 minute.** Discard the flow-through.
10. **Dry silica membrane.** Centrifuge at 14.000 rpm for 3 minutes.
11. Place the MicroSpin Column into a 1.5 mL nuclease-free tube (not provided) and add **50-100 µl Pre-heat the Elution Buffer** at 70°C on the centre of the white filter membrane. Incubate **at room temperature for 2 minutes.**
12. **Centrifuge** the spin column-tube assembly **at 14.000 rpm for 1 minute,** then discard the column. The DNA is now ready for downstream applications.

#### **4. PROBLEM GUIDE AND POSSIBLE ANSWER**

For any doubts or additional questions about the protocol, please contact the technical service of DANAGEN-BIOTED S.L [info@danagen.es](mailto:info@danagen.es)